



no to NECK PAIN

...no to TEXT-NECK

Ghayth A., Leah B., Lauren M., Veronica M.

01

BACKGROUND

Cervical flexion causes non-specific neck pain.

03

DATA FLOW

Building prototype and assessing functionality.

02


SOLUTION

Design proposition: features & constraints.

04

LAWS & ETHICS


Implementation plan for final product.



01

BACKGROUND

How does **cervical flexion** lead to non-specific neck pain?



NON-SPECIFIC NECK PAIN

A global public health burden.



UNDEFINED:

no generally accepted medical definition; no root cause.

PREVALENT:

afflicts ~**71%** of the population at some point in their lives.

DESTRUCTIVE:

associated with loss of ROM and decreased QoL.

TEXT-NECK

→ Dubbed as 'text-neck' only recently, **cervical spine flexion** is a common posture for individuals using *computers, tablets* and *smartphones*.

→ This posture has since been flagged as a **major risk factor** for developing **nonspecific neck pain** due to *increased loading*.



To prevent the development of nonspecific neck pain resulting from cervical spine flexion, a **monitor and alarm system** is needed to detect and correct 'text-neck' posture.

—**NEED STATEMENT**



DESIGN CRITERIA

MUST BE: **constraints**

SHOULD BE: **features**

CONSIDER: **options**




Accurate	Algorithm-Efficient	Sensor placement/DAQ
Wearable & Safe	Easy to conceal	Hardware, housing needs
Communicative	Durable	Buzzer type/placement
Reproducible	Easy to use	Sensor placement

02

SOLUTION

A sensor-equipped **balaclava** to monitor neck posture and alert user if sustained flexion is detected



<50 deg: cervical flexion

[9]

MONITORING SYSTEM

angle between C7 and the ear relative to horizontal



SENSOR →
AXIS OF EAR

Captures **signal** (end effector of cervical spine)

SENSOR →
C7 SPINE


Provides **reference** (eliminating noise)

NOISE → **L-T SPINAL LEVELS**

Eliminated by placing reference at C7

NOISE → **HIP & BELOW**

Eliminated by placing reference at C7



BALACLAVA DESIGN

REMOVEABLE

Cloth at front can come on and off for optimal user comfort

HOUSING

Contains (1) breadboard with **micro** circuit, (2) **power** supply and (3) **LEDs**; hangs from straps



BUZZER

Placed in the breadboard; close enough to hear


SENSOR 1

BNO055 signal at axis of ear


SENSOR 2

BNO055 reference at C7


MATERIALS



Breadboard



ESP 32



BNO550




Buzzer



LED



M/F Jumper
Wire



QWIC Wire




USB to
micro-USB




Balaclava

HARDWARE SETUP




HOUSING DESIGN




03

DATA FLOW


How data flows through our system... from acquisition to visualization



SYSTEM BLOCK DIAGRAM



DATA FLOW (DAQ)



OBTAINING DEGREE OF FLEXION

```
76 float angle=(event.orientation.z);//angle of the ear relative to the horizontal
77 bno2.getEvent(&event2);
78 float reference=((event2.orientation.z)+50);// angle of C7 relative to the horizontal
79
80 if (i<20)
81 {
82     anglesum=anglesum+angle;
83     refsum=refsum+reference;
84 }
85 else {
86     float angleavg=anglesum/20;
87     float refavg=refsum/20;
88     float sum = (refavg - angleavg);
89     SerialBT.print((90-sum));
```

1

Calculate **average reference** angle and **ear** angle


2

Calculate **neck** angle

3

Display neck angle **relative to the horizontal**


DATA FLOW (STATE MACHINE)



DATA VISUALIZATION

	A	B
1	Average angle	25.93777399
2	Time	Neck Angle
3	1.04	80.43
4	1.04	78.79
5	1.29	76.67
6	1.29	74.75
7	1.55	72.65
8	1.55	70.66
9	1.8	68.55
10	1.8	66.4
11	2.05	63.95
12	2.05	61.45
13	2.29	58.54

EXCEL FILE




PLOT

04

USER IMPACT

A review and **assessment** of our
device build & user feedback
against initial criteria



PATIENT FEEDBACK

Patient neck angle >50 degrees

Buzzer OFF
Green LED ON

Patient neck angle <50 degrees for
 <30 seconds **before** Red LED
OR

Patient neck angle >50 degrees for
 <30 seconds **after** Red LED

Turn buzzer
OFF; Yellow
LED ON


Patient neck angle <50 degrees for
 >30 seconds

Buzzer ON
Red LED ON

CLINICIAN FEEDBACK

OUTPUTS AFTER WEARABLE SESSION:

- Neck angle (red)
 - ◆ Gives visual feedback on neck angle
- User compliance (green)
 - ◆ Good compliance = 0 on the plot
 - ◆ Bad compliance = peaks on the plot
 - ◆ Compliance = how quickly the user responds to live feedback
 - ◆ The shorter the compliance peak the faster the response to poor posture
- Plotted against the reference threshold = 50 degrees (blue)




CALIBRATION & VALIDATION

Use the **built in calibration** of the gyroscope



Calibrate so a normal cervical neck position is 90 degrees relative to the horizontal in the sagittal plane



Use tests to **validate accuracy** of design in both a **simulation and clinically**

Set both the reference and angle sensors to 90 degrees relative to the horizontal

Monitor angle measurements and relative angle (neck angle)

Move the angle sensor to various desired and known positions

Monitor angle measurements and relative angle (neck angle)

HOW DESIGN MEETS CRITERIA

MUST BE: **constraints**

SHOULD BE: **features**

CONSIDER: **options**



Accuracy **validated**

Simple/**efficient**
code written

Code outputs correct
neck angle measures



Balaclava

Small components

Secure **housing** designed



Buzzer/LED **alert**
system

Durable/standard
equipment

Buzzer/LEDs placed in
housing



Reproducible

Simple **wearable**

In line with the **ear** & **C7**


05

USER/ETHICAL CONSIDERATIONS



CLASS 1 MEDICAL DEVICE

- Our balaclava aims to monitor neck posture to combat 'text-neck' and prevent the onset of non-specific neck pain.
- Because its **intended use is to prevent a specific disease**, i.e. non-specific neck pain, this device would be considered, under both Health Canada and the FDA, a **medical device**.
- It would be considered a **Class I** device because it is very low-risk and simple in design.



CAD/US. CLASS 1 REGULATIONS

FDA: Gen. Controls Only

requires adherence to manufacturing (**GMP**) and **labelling** standards



HEALTH CAD: Premarket Approval (PMA)

medical Device License (**MDL**) which has an associated cost of **C\$389**, per application

LIABILITY & PATENT PROTECTION

Device provides **inaccurate** data

We engineers are held **liable**

Thus we conduct **testing** before market release

Device is **similar** to others on the market

Thus we conduct a search of existing **patents**

And we build upon our design **uniquities**

Alters clinical decisions on the basis of falsehoods (**dangerous**)

Device leaks **private** data without consent or due to **poor security**

We engineers are held **liable**

Thus we ensure **security** and obtain **explicit consent**

Transparent about data collection

EQUITY, INCLUSIVITY & DIVERSITY



Belonging

Feeling accepted as member of **society**, valued, connected

Inclusivity

Implement design features to accommodate **disabilities**; e.g., for patients with reduced motor function



Equity

Ensure our device does not perpetuate social inequality (e.g. via insurance companies) and **stigmatisation**



Diversity


Ensure our device is **accessible** and beneficial to a diverse group of users, outreach extending to **marginalised groups** (emphasize **affordability**)



06

FINAL DEMO

A live demonstration of our
device in action



DEMO (VIDEO)

3p04_final.mov



BILL OF MATERIALS - COSTS

Item	Cost	QTY	Sourced From	Total
Breadboard	\$4.45	1	Digikey	\$4.45
ESP 32	\$15.34	1	Digikey	\$15.34
BNO055	\$45.95	2	Digikey	\$91.90
Buzzer	\$2.92	1	Digikey	\$2.92
LED	\$0.33	3	Digikey	\$1.00
330 Ω resistor	\$0.10	3	Digikey	\$0.30
USB2.0 to Micro-USB Cable	\$7.16	1	Digikey	\$7.16
M/F 6" Jumper Wire	\$3.03	1	Digikey	\$3.03
Jumper Wire M to M 3" 28AWG	\$1.95	1	Digikey	\$1.95
QWIIC Wires	\$2.16	2	Digikey	\$4.32
Micro-USB Wireless Charger	\$7.59	1	Amazon	\$7.59
Balaclava	\$11.74	1	Amazon	\$11.74
3D Printed Housing	\$0.77	1	N/A	\$0.77
TOTAL				\$152.47

ACQUISITION COSTS

DIGIKEY ORDER

Before tax: \$129.12

After tax: \$145.90

Shipping free in
Canada for >\$100

\$145.90

AMAZON ORDER

Before tax: \$19.33

After tax: \$21.84

Shipping & customs
fee = \$14.99

\$36.83

REFERENCES

- [1] R. Fejer, K. O. Kyvik, and J. Hartvigsen, "The prevalence of neck pain in the world population: A systematic critical review of the literature," *European Spine Journal*, vol. 15, no. 6, pp. 834–848, 2005.
- [2] "A standardized approach to the diagnosis and treatment of Spinal Pain," *Low Back and Neck Pain*, pp. 190–226, 2004.
- [3] N. F. Mahmoud, K. A. Hassan, S. F. Abdelmajeed, I. M. Moustafa, and A. G. Silva, "The relationship between forward head posture and Neck Pain: A systematic review and meta-analysis," *Current Reviews in Musculoskeletal Medicine*, vol. 12, no. 4, pp. 562–577, 2019.
- [4] L. Hofstetter, M. Häusler, P. Schweinhardt, U. Heggli, D. Bron, and J. Swanenburg, "Influence of axial load and a 45-degree flexion head position on cervical spinal stiffness in healthy young adults," *Frontiers in Physiology*, vol. 12, 2021.
- [5] A. N. Vasavada, D. D. Nevins, S. M. Monda, E. Hughes, and D. C. Lin, "Gravitational demand on the neck musculature during tablet computer use," *Ergonomics*, vol. 58, no. 6, pp. 990–1004, 2015.
- [6] W. Y. Wong and M. S. Wong, "Trunk posture monitoring with inertial sensors," *European Spine Journal*, vol. 17, no. 5, pp. 743–753, 2008.
- [7] D. Giansanti, M. Dozza, L. Chiari, G. Maccioni, and A. Cappello, "Energetic assessment of trunk postural modifications induced by a wearable audio-biofeedback system," *Medical Engineering & Physics*, vol. 31, no. 1, pp. 48–54, 2009.
- [8] L. Simpson, M. M. Maharaj, and R. J. Mobbs, "The role of wearables in Spinal Posture Analysis: A systematic review," *BMC Musculoskeletal Disorders*, vol. 20, no. 1, 2019.
- [9] K. Dunleavy, J. Neil, A. Tallon, and D. E. Adamo, "Reliability and validity of cervical position measurements in individuals with and without chronic neck pain," *The Journal of manual & manipulative therapy*, vol. 23, no. 4, pp. 188–196, 2015.